

2010

HEAD TO TOE PROTECTION WITH

DermaSaver™




Revolutionary Skin Protection & Pressure Reduction

Unlike foams, gels, fleece, and polyfil, DermaSaver™ reduces ALL of the conditions that cause skin breakdown.




LEFT SHIN 10/1/07



LEFT SHIN 10/4/07



LEFT SHIN 10/11/07



LEFT SHIN 10/16/07



“ All the DermaSaver products are soft & springy. I’ve seen a great reduction of skin tears on my unit. These products are great for prevention from breakdown for heel, knees & elbows. ”
~ Nursing Staff, Greenbriar Terrace, Nashua, NH

CALL HIPSAVER AT 1-800-358-4477



**DermaSaver
by HipSaver Inc.**
7 Hubbard Street
Canton, MA 02021
www.dermasaver.com
email: cs@dermasaver.com

TABLE OF CONTENTS

4	DermaSaver SkinTubes for Legs
7	DermaSaver SkinTubes for Arms
10	DermaSaver Heel & Ankle Protection
12	DermaSaver Pressure Reduction & Protection
15	HipSaver Hip Protectors

Dear Health Care Provider,

We developed DermaSaver three years ago to provide a superior preventative product for skin tears, bruising, and pressure ulcers. After listening to care-givers telling us what they wanted and what was missing from the products they had available before DermaSaver, we built in the following features:

- ◆ Products that fit right and stay positioned on the desired area.
- ◆ Smooth fine fiber knit surface yarns to reduce friction and shear
- ◆ Soft & cushy product with the same resillance of skin (Indentation load deflection)
- ◆ Antimicrobial properties
- ◆ Does not adhere to wounds
- ◆ Breathable fabric that does not trap moisture
- ◆ Totally launderable

To our surprise, these features not only created a preventative effect but also a therapeutic effect on existing wounds. All of these features in combination allow existing wounds to heal by protecting from additional environmental assaults.

Throughout this catalog you will notice before and after photos from our Case Study Collection. You can view all the case studies in full detail at www.dermasaver.com or ask us to mail you our 20 page Case Study Collection Brochure.

Please take the time to review these impressive case studies and their implications to Improve your care plan. The prevention of a single would improves life quality and saves real dollars for your facility. And surveyors have been impressed with the successful results achieved with DermaSaver.

As always, please call me with questions, comments or ideas.

Sincerely,

Take a look at the **Stay-Put™ Heel Protectors** on Page 10. They're a hot seller!




PS. We're always striving to provide something new and now all DermaSaver products are available with a pure silver thread lining to promote healing, reduce pain and kill bacteria. See page 3 for details.

Ed Goodwin (right) Founder and President of HipSaver Company and inventor of DermaSaver products enjoys a Spring day with his rescued Rottweiler Sonja. Check out Sonja's favorite website 4theloveofdog.org Sonja thanks you for supporting your local animal rescue organizations.

The unique, patent pending constructions of DermaSaver products were developed by Ed Goodwin, founder and inventor of HipSaver hip protectors. DermaSaver products are manufactured exclusively by HipSaver, Inc.

Call 1-800-358-4477

Get Head to Toe Protection with **DermaSaver™**



REVOLUTIONARY MICROSPRING™ TEXTILE OFFERS BREATHABILITY, MOISTURE TRANSFER AND ANTI-MICROBIAL PROTECTION NOT POSSIBLE WITH GELS, FOAMS OR FLEECE.


Introducing the DermaSaver Silver Line

All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. The X-Static fiber has a layer of pure silver permanently bonded to it. It is all natural, non-toxic and allows for heat transfer. Silver is the most conductive metal and because of its conductivity, it promotes healing and relieves pain. Silver has of course been used in other medical applications for years: silver nitrate in infant's eyes, Silverdene cream for burns and Colloidal Silver. And perhaps most importantly, the X-Static fiber has been shown to kill 99.9% of bacteria in less than one hour of exposure. It is permanent and will not wash out.

WHAT MAKES DERMASAVER WORK SO WELL? THE SCIENCE BEHIND THE RESULT.



Durometer reading from skin (30)



The Aegis Microbe Shield inhibits microorganism growth on treated area ONLY. There is no microbial encroachment and no leaching outside the treated area.



Durometer reading from DermaSaver (30) matches skin

Microfilaments flex to reduce pressure on bony prominences such as the heel, elbow, vertebra and trochanter and to protect arms and legs from bumps and other traumas that cause skin tears and bruises. DermaSaver SkinTubes are made of two layers of MicroSpring Textile using 8000 microfilaments per square inch to act as mini pressure reducing springs. DermaSaver pressure reduction products have up to 20 thousand microfilaments per square inch.

DermaSaver MicroSpring Textile is designed to have the same "cushiness" (durometer compression resistance) as human skin. It is like adding an external layer of skin and subcutaneous fat. The surface yarns are a fine knit which reduces shear and friction.

DermaSaver MicroSpring Textile is 94% Polyester and 6% Spandex, Latex-free and meets the CA-117 flame retardant standard.

All layers of DermaSaver products are protected by the Aegis Microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product only. It is permanently bonded to the fibers, does not leach, or volatilize. It does not foster adaptive micro-organisms.



DermaSaver SkinTubes™ for Legs

- ◆ SUPERIOR PROTECTION FROM SKIN TEARS
- ◆ NON-CONSTRICTING
- ◆ ANTI-MICROBIAL
- ◆ BREATHABLE
- ◆ TOTALLY LAUNDERABLE

Two layers of soft, lightweight MicroSpring Textile™ form a loose-fitting tube that stays in place to protect fragile skin from pressure, friction and minor traumas that result in tears, bruising and minor abrasions. Unlike traditional fabric stockings, which are tight and binding, the layers of microfilaments allow air-flow and moisture transfer.

Simply slip the tube over the leg, there are no hook and loop fasteners or tight fitting spandex. The narrower end of the tube gently rests on the ankle to keep it from falling off. The cross microfilaments in the MicroSpring Textile allows these soft tubes to be self-supporting – they won't slip down or bunch up.

SkinTubes are protected by the Aegis Microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product.

You can use creams and ointments with SkinTubes. All SkinTubes can be completely washed and dried even in high temperature institution laundries.

SkinTubes are also available with double protection for knees.



THESE SOFT FLEXIBLE TUBES FIT LOOSELY AND REMAIN ELONGATED ON ARMS & LEGS

TWO LAYERS OF MICROSPRING TEXTILE PROTECT FRAGILE SKIN WITH 8000 MICROFILAMENTS PER SQUARE INCH.



Full Leg Tube

Two layers of MicroSpring Textile protect from the ankle to above the knee.

- LT1300XS** (X-Small)
- LT1300S** (Small)
- LT1300M** (Medium)
- LT1300L** (Large)
- LT1300XL** (X-Large)

Full Leg Tube with Double Knee

Four layers of MicroSpring Textile protect the knee. Recommended for non-ambulatory residents only.

- EKT1300XS** (X-Small)
- EKT1300S** (Small)
- EKT1300M** (Medium)
- EKT1300L** (Large)
- EKT1300XL** (X-Large)



Knee Tube

Two layers of MicroSpring Textile protects the knee.

- KT1350XS** (X-Small)
- KT1350S** (Small)
- KT1350M** (Medium)
- KT1350L** (Large)
- KT1350XL** (X-Large)

Double Knee Tube

Four layers of MicroSpring Textile protects the knee.

- 4KT1350XS** (X-Small)
- 4KT1350S** (Small)
- 4KT1350M** (Medium)
- 4KT1350L** (Large)
- 4KT1350XL** (X-Large)

SIZING FOR LEG TUBES:

Measure around the widest part of the calf. Choose the larger size if measurement is between sizes:

- X-Small:** 7"- 9"
- Small:** 9"- 12"
- Med:** 12"- 15"
- Large:** 15"- 19"
- X-Large:** 19"- 23"

All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).

CASE STUDY SHOLOM HOME WEST, ST. LOUIS PARK, MN

BEFORE USING DERMASAVER SHIN-KNEE TUBES



10/11/07



Shin~Knee Tube

Designed especially for residents in wheelchairs, the Shin~Knee Tube protects both the shin & the knee with 2 layers of MicroSpring textile. The pleated contoured cup stays in position on the knee.

- SKT1350XS** (X-Small)
 - SKT1350S** (Small)
 - SKT1350M** (Medium)
 - SKT1350L** (Large)
 - SKT1350XL** (X-Large)
- same sizing as leg tubes

AFTER 21 DAYS USING DERMASAVER SHIN-KNEE TUBES



11/1/07

We had one resident with Peripheral Vascular Disease (PVD) who was getting bruises that would bleed under the skin and take the skin off her shins. Since using the DermaSaver Shin Tubes, her shins have healed and she has not had another bruise or tear.

~ Dee Hendricks, DON, HillCrest Healthcare Service, Hawarden, IA



Shin Tube

Two layers of MicroSpring Textile protect from the ankle to below the knee

STI350XS (X-Small)

STI350S (Small)

STI350M (Medium)

STI350L (Large)

STI350XL (X-Large)

same sizing as leg tubes

Amputee "Stump" Cover

A contoured cup of 2 layers of MicroSpring Textile provides protection for the amputation site. Stretch fabric helps to keep the cover in place.

SCI350XS (X-Small)

SCI350S (Small)

SCI350M (Medium)

SCI350L (Large)

SCI350XL (X-Large)

The DermaSaver arm & leg tubes are great products for skin protection. Residents find them soft and comfortable, and not tight like other products we've tried. They work really well with residents who are self mobile in the wheelchair and have a tendency to bump their arms and legs. We really like your products.

~ **Loretta Gerving, DON, Marian Manor HealthCare Center, Glen Ulin, ND**

SIZING FOR AMPUTEE STUMP COVER:

Measure the circumference 2" above the amputation site. Choose the larger size if measurement is between sizes:

X-Small: 8"-10"

Small: 10"-12"

Med: 12"-15"

Large: 15"-18"

X-Large: 18"-22"

CASE STUDY

WESBURY UNITED METHODIST COMMUNITY



All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).

DermaSaver SkinTubes™ for Arms



- ◆ SUPERIOR PROTECTION FROM SKIN TEARS
- ◆ NON-CONSTRICTING
- ◆ ANTI-MICROBIAL
- ◆ BREATHABLE
- ◆ TOTALLY LAUNDERABLE

Forearm Tube

Two layers of MicroSpring Textile protect from the wrist to below the elbow.

- FTI350XS** (X-Small)
- FTI350S** (Small)
- FTI350M** (Medium)
- FTI350L** (Large)
- FTI350XL** (X-Large)

Two layers of soft, lightweight MicroSpring Textile™ form a loose-fitting tube that stays in place to protect fragile skin from pressure, friction and minor traumas that result in tears, bruising and minor abrasions. Unlike traditional fabric sleeves, which are tight and binding, the layers of microfilaments allow air-flow and moisture transfer.

Simply slip the tube over the arm, there are no hook and loop fasteners or tight fitting spandex. The narrower end of the tube gently rests on the wrist to keep it from falling off. The cross microfilaments in the MicroSpring Textile allows these soft tubes to be self-supporting - they won't slip down or bunch up.

SkinTubes are protected by the Aegis Microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product.

You can use creams and ointments with SkinTubes. All SkinTubes can be completely washed and dried even in high temperature institution laundries. SkinTubes are also available with double protection for elbows.



Arm Tube

Two layers of MicroSpring Textile protects from the wrist to above the elbow.

- ATI300XS** (X-Small)
- ATI300S** (Small)
- ATI300M** (Medium)
- ATI300L** (Large)
- ATI300XL** (X-Large)

Full Arm Tube with Double Elbow

Full arm protection with four layers of MicroSpring Textile covering the elbow. Replaces traditional foam elbow protectors.

- EET1300XS** (X-Small)
- EET1300S** (Small)
- EET1300M** (Medium)
- EET1300L** (Large)
- EET1300XL** (X-Large)

CASE STUDY RIDGEWOOD CARE CENTER, RACINE, WI



THESE SOFT FLEXIBLE TUBES FIT LOOSELY AND REMAIN ELONGATED ON ARMS & LEGS



Forearm Tube With Knuckle Protector

Two layers of MicroSpring Textile provide all-in-one protection from the knuckles to below the elbow.

- FTKPI3XS** (X-Small)
- FTKPI3S** (Small)
- FTKPI3M** (Medium)
- FTKPI3L** (Large)
- FTKPI3XL** (X-Large)



Full Arm Tube With Knuckle Protector

Two layers of MicroSpring Textile provide all-in-one protection from the knuckles to above the elbow.

- ATKPI3XS** (X-Small)
- ATKPI3S** (Small)
- ATKPI3M** (Medium)
- ATKPI3L** (Large)
- ATKPI3XL** (X-Large)

I LOVE YOUR PRODUCT. MOM HAS HAD TERRIBLE SKIN TEARS FOR THE LAST COUPLE OF YEARS AND THIS IS THE FIRST PRODUCT THAT HAS ACTUALLY PREVENTED THEM! BLESS YOU ALL.

~ Kit Frazer, Altamonte Springs, FL



→ Arm~Bow Tube

Two layers of MicroSpring textile protect the forearm while 5 layers are kept in place to protect the elbow. Will not roll around the elbow.

- ABTI350XS** (X-Small) *same sizing as arm tubes*
- ABTI350S** (Small)
- ABTI350M** (Medium)
- ABTI350L** (Large)
- ABTI350XL** (X-Large)



ZipTube

Two layers of MicroSpring Textile protects the arm. Zips open to allow for inspection of IV site and wound location. Loose enough to rotate. Zipper doesn't come in contact with skin.

- ZTI300-XS** (X-Small)
- ZTI300-S** (Small)
- ZTI300-M** (Medium)
- ZTI300-L** (Large)
- ZTI300-XL** (X-Large)



SIZING FOR ARM TUBES:

Measure around the widest part of the forearm. Choose the larger size if measurement is between sizes:

- X-small:** 5"- 7"
- Small:** 7"- 10"
- Med:** 10"- 12"
- Large:** 12"- 15"
- X-large:** 15"- 18"

All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).

CASE STUDY

NATHAN LITTAUER HOSPITAL
NURSING HOME, GLOVERSVILLE, NY

BEFORE USING DERMASAVER ELBOW TUBE



OCT. 03, 2007: 88 year old female with skin tear above right elbow. “V” shaped 1.25cm x 0.25cm scant amount of bleeding. Dressed with Mepitel (Molnlycke), Bacitracin, DSD.

AFTER 14 DAYS USING DERMASAVER ELBOW TUBE



OCT. 16, 2007: Healed, closed wound. Continued use of Remedy Skin Repair Cream and DermaSaver Elbow Tube for prevention.



Elbow Tube

Two layers of MicroSpring Textile protects the elbow

ET1350XS (X-Small)

ET1350S (Small)

ET1350M (Medium)

ET1350L (Large)

ET1350XL (X-Large)

Double Elbow Tube

Four layers of MicroSpring Textile cover the elbow.

Replaces traditional foam elbow protectors.

4ET1350XS (X-Small)

4ET1350S (Small)

4ET1350M (Medium)

4ET1350L (Large)

4ET1350XL (X-Large)

“We’ve had great success with your DermaSaver products. On one resident in particular, the DermaSaver arm tubes not only prevented skin tears but healed her existing tears much quicker than anything we had tried before.”

~ Diane Husar, Dir. Quality Assurance, United Methodist Homes, Binghamton, NY

DermaSaver Heel and Ankle Protection



Stay-Put™ Heel Protector

Five layers of MicroSpring Textile cover the heel to reduce pressure and friction on the heel and malleolus. Unique design keeps the heel protector from rolling around or slipping off the foot. Comfortable enough to wear all day and all night. *Can be worn while walking.*

HPI300-S (Small) **HPI300-M** (Medium)
HPI300-L (Large) **HPI300-XL** (X-Large)

See sizing diagram on page 11



Stay-Put™ Plus Heel Protector with toe covering

The same Stay-Put Heel Protector as above with fold-over flap to cover the toes. Same sizes as Stay-Put Heel Protector. Available in regular and long lengths.

HPPI300-S (Small) **HPPI300-M** (Medium)
HPPI300-L (Large) **HPPI300-XL** (X-Large)

See sizing diagram on page 11



Relevator™ Heel Elevator

Two inches of contoured, cushy polyfil covered by MicroSpring Textile lift the heel in total comfort. Full circle coverage provides consistent elevation during movement. One-piece, slip-on with no fasteners, also covers both malleoli. Fits loosely yet will not fall off. Durable and institution launderable for economical repeated long-term and non patient-specific use. Replaces bed pillow elevation.

REL300-S/M (Small/Medium) **REL300-L/XL** (Large/X-Large)

See sizing diagram on page 11



Because it is multi-layered, the Relevator may take over one hour to fully dry depending upon dryer temperature. This product will withstand the institution high drying temperatures. The Relevator Heel Elevator is designed to be worn by persons while laying down. Persons who raise their knees and have their foot flat on the bed will not benefit from the Relevator. Check each individual for proper fit and applicability of this product. If you have any questions call Ed Goodwin at 1-800-358-4477.

CASE STUDY

SENIOR VILLAGE PERRYTON, TX

Launder Guard

BOTH THE DERMASAVER STAY-PUT AND STAY-PUT PLUS HEEL PROTECTORS COME WITH THE LAUNDER GUARD. SIMPLY CLOSE THE SNAP PRIOR TO LAUNDERING TO KEEP THE VELCRO CLOSURES FROM SNAGGING GARMENTS IN THE LAUNDRY.

BEFORE USING DERMASAVER STAY-PUT HEEL PROTECTOR



AFTER USING DERMASAVER STAY-PUT HEEL PROTECTOR



All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).

Sizing for Heel Protection:

To determine the correct size, measure around foot & heel as indicated in the diagram.

Stay-Put™ Heel Protector

Small: Less than 12"

Medium: 12" – 14 1/2"

Large: 14 1/2" – 15 1/2"

X-Large: Greater than 15 1/2"

Relevator™ Heel Elevator

Small/Medium: Less than 13"

Large/X-Large: Greater than 13"



Measure the widest distance all the way around the heel/instep



Ankle Bumpers™

Three layers of MicroSpring Textile wrap around the ankle to protect the malleoli.

Measure around widest part of ankle bone.

ANK13-S (Small) 9" or less

ANK13-M (Medium) 9"-11"

ANK13-L (Large) 11"-13"

ANK13-XL (X-Large) 13" or greater

CASE STUDY

FLORIDA HOSPITAL WOUND CLINIC, ORLANDO, FL

BEFORE USING DERMA SAVER RELEVATOR



12/20/07

AFTER 15 DAYS USING DERMA SAVER RELEVATOR



1/4/08

Within 4 days of applying the Stay-Put Heel Protectors, the redness from a stage 1 pressure ulcer was resolved on a bedfast resident with upper and lower body contractures.

~ Nursing Staff, Greenbriar Terrace, Nashua, NH

DermaSaver Pressure Reduction & Protection

Layers of MicroSpring Textile put thousands of microfilaments per square inch to act as mini pressure-reducing springs. Microfilaments flex to reduce pressure on bony prominences such as the heel, elbow, vertebra and trochanter. Ample ventilation space between the microfilaments and textile layers allow for air flow and moisture transfer.

DermaSaver products are protected by the Aegis Microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product.

You can use creams and ointments with DermaSavers.

All DermaSaver products can be completely washed and dried even in high temperature institution laundries.



WHEELCHAIR TOP CUSHION HAS FABRIC CORNER TIES SO THE PAD CAN BE USED ON THE SEAT OR WHEELCHAIR BACK. TIES CAN BE CUT OFF IF NOT NEEDED.

← Pressure Reduction Interface Pads

Use over mattress, wheelchair cushion or other surfaces. Five layers of MicroSpring Textile provides pressure reduction and protection from friction and shear. Perfect for torso areas such as scapula, sacrum, ribs, trochanter, ischial tuberosity, condyles and heels.

FP1818 Wheelchair Top Cushion (18" X 18")

FP3524 Area Flat Pad (35" X 24")

FP3578 Mattress Topper (35" X 78")

We've had problems with other protectors sticking to the wounds, but not with DermaSaver. I really like the air permeability and antimicrobial qualities of the products, too.

~ Diane Husar, Dir. Quality Assurance, United Methodist Homes – Elizabeth Church Campus Binghamton, NY

CASE STUDY

Citronelle Convalescent Center, Citronelle, AL

BEFORE USING DERMA SAVER DERMA BRIEF



12/10/07 Sacral Wound – STAGE II

1 WEEK USING DERMA SAVER DERMA BRIEF



12/17/07 STAGE II

2 WEEKS USING DERMA SAVER DERMA BRIEF



12/24/07 STAGE II HEALING

AFTER 4 WEEKS NO OPEN AREAS TO SACRAL WOUND



1/4/08 Wound HEALED

All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).



DermaBrief™

Three layers of MicroSpring Textile cover the Sacral area to reduce pressure and shear. Also available covering the hips and sacral area. Sizes based on hip circumference Choose Pull-Up or Snap Model.

Pull-Up Model, Hip & Sacral area:

- DBI3-XS** (28"-31")
- DBI3-S** (32"-35")
- DBI3-M** (36"-39")
- DBI3-L** (40"-44")
- DBI3-XL** (45"-50")
- DBI3-XXL** (51"-57")

Snap Model, Hip & Sacral area:

- DBWS-XS** (28"-31")
- DBWS-S** (32"-35")
- DBWS-M** (36"-39")
- DBWS-L** (40"-44")
- DBWS-XL** (45"-50")
- DBWS-XXL** (51"-57")

FRONT VIEW OF SNAP MODEL



Pull-Up Model, Sacral area only:

- DBTI3-XS** (28"-31")
- DBTI3-S** (32"-35")
- DBTI3-M** (36"-39")
- DBTI3-L** (40"-44")
- DBTI3-XL** (45"-50")
- DBTI3-XXL** (51"-57")

Snap Model, Sacral area only:

- DBWST-XS** (28"-31")
- DBWST-S** (32"-35")
- DBWST-M** (36"-39")
- DBWST-L** (40"-44")
- DBWST-XL** (45"-50")
- DBWST-XXL** (51"-57")



DermaWrap™

Single layer of MicroSpring Textile (5' x 10') provides flexible and even compression wrap for swollen limbs. End can be secured by self-gripping tape or other means. Can be cut to desired length.

DWI300
Per Roll (5'x10')



Oxygen Tube Covers

One layer of MicroSpring Textile with pressure-sensitive medical grade adhesive covers oxygen tubes to reduce friction against skin, especially between the back of the ear.

FTCI300 Package of 10 Covers (2³/₄" x 1¹/₂")

Note: Because of the adhesive used on them, these are not launderable.



Pressure Reduction Pillow Case

Standard pillowcase with 5 layers of MicroSpring Textile on top side to reduce pressure and shear on back of head and ears while sleeping or laying down.

PCI300 Standard Size (26" x 20")



Derma-Tee™

Tee-shirt has a full back panel of 2 layers of MicroSpring Textile. Reduces shear and pressure on scapula and vertebral bones. Excellent for those with kyphosis.

Sizes based on chest circumference:

- DTI3S** (34"-36")
- DTI3M** (38"-40")
- DTI3L** (42"-44")
- DTI3XL** (46"-48")
- DTI3XXL** (50"-52")

We are very impressed with the DermaSaver products. The resident that used the Pressure Reduction pillowcase commented that she had not slept so well in a long time. We put the Stay-Put Heel Protector on a gentleman that would never keep anything on but he has been wearing this product consistently. Your products are very nice!

~ Diane Kiel, Rehab Coord
Cambridge Health Center, Cambridge, MN



Finger Separator

Four soft layers of MicroSpring Textile are placed between each finger and the cushioning palm pillow provides maximum comfort. Adjustable soft band with Velcro®.

Will not flatten out.

FSI300 One Size



Palm Pillow

Layers of MicroSpring Textile rolled to 1 1/2" thick pillow prevents digging fingernails into palms. Adjustable soft band with Velcro®.

PPI300-S (Small)

PPI300-L (Large)



DermaMitt™

Two layers of MicroSpring Textile with an inner core layer of half-inch polyfil helps prevent scratching and tube pulling. Protects fragile skin on hands. Velcro® securing band is folded under soft cuff.

HMI300-SM (Sm-Med)

HMI300-LX (Lrg-XL)



Knuckle Protector

Two layers of MicroSpring Textile cover sensitive knuckles to prevent bumping and scraping. Perfect for those using wheelchairs.

KPI300SM (Small-Medium)

KPI300LX (Large - X-Large)



CASE STUDY

Brandywine Nursing and Rehab Center, Wilmington, DE

BEFORE USING DERMASAVER KNUCKLE PROTECTOR



1/11/08

AFTER USING DERMASAVER KNUCKLE PROTECTOR



2/4/08



Long Knuckle Protector

Protects more of the finger area.

KPLI300SM (Small-Medium)

KPLI300LX (Large-X-Large)

All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).

REDUCE FALL RELATED HIP FRACTURES

HipSaver Hip Protectors

Protects vulnerable hips during a fall.

- ◆ SOFT, SEWN-IN, TAPERED PADS
- ◆ PROVEN IN BIOMECHANICAL TESTING AND CLINICAL STUDIES
- ◆ MACHINE WASH & DRY EVEN IN HOT INSTITUTION LAUNDRIES

HipSaver Nursing Home Model for easy dressing

Wide leg holes, roomy fit, easy to pull up and down.



Nursing Home



Nursing Home with male fly-front



Add Tailbone pad for

HipSaver QuickChange™ and HipSaver Wrap&Snap for incontinence management

QuickChange has a drop down panel with only 2 snaps for easy changing of incontinence products. Wrap&Snap snaps completely off to lay flat for hard to dress residents



QuickChange



Wrap&Snap



Add Tailbone pad for

HipSaver SlimFit™ for outpatients and fashion conscious residents

Sleek design – discreet under clothing. Ideal for fashion conscious residents and outpatients



SlimFit



SlimFit with male fly-front



Add Tailbone pad for

All HipSavers come in 6 sizes, X-Small thru XX-large. Measure around the widest part of the hip to determine correct size.



HipSaver OpenBottom™

Open crotch design for easy self-toileting. Use alone or residents can wear their own underwear over it.

HipSaver Open-Bottom 3-Snap™

Designed to be worn over incontinence products. Open crotch area reduces chance of soiling the HipSaver resulting in less frequent changing.



HipSaver SoftSweats™

Proven HipSaver protection in soft comfortable sweatpants and shorts.

- Six hip sizes and 2 lengths
- No elastic or zippers at ankle
- Choose navy or grey

Available with:

- ¹ Hip pads:
- ¹ Hip & Tailbone pad:
- ¹ Hip & Knee pads:
- ¹ Hip, Tailbone & Knee pads:



All HipSavers are available with EZ Pull-Up Handles to aid those with weakened grip strength.


2010

DermaSaver™

DERMASAVER PRESSURE REDUCTION & PROTECTION



All DermaSaver products are now available with X-Static® silver thread knit into the layer closest to the skin. Add 'SIL' after product #. (See page 3 for details).



Revolutionary MicroSpring™ Textile Offers Breathability, Moisture Transfer and Anti-Microbial Protection not possible with gels, foams, fleece or polyfil.

UNLIKE FOAMS, GELS, FLEECE AND POLYFIL, DERMASAVER™ REDUCES ALL OF THE CONDITIONS THAT CAUSE SKIN BREAKDOWN.

Breathability and Moisture Transfer

Ample ventilation space between the microfilaments and textile layers allow for air flow and moisture transfer.

MicroSpring Textile

Microfilaments flex to reduce pressure on bony prominences such as the heel, elbow, vertebra and trochanter. Layers of MicroSpring Textile put thousands of microfilaments per square inch to act as mini pressure reducing springs.

Protected by the ÆGIS Microbe Shield

All layers of DermaSaver products are protected by the Aegis Microbe shield which uses an EPA registered antimicrobial technology to inhibit the growth of micro-organisms on the product. It is permanently bonded to the fibers, does not leach, or volatilize.

Totally Launderable

Machine wash & dry even in the high temperature of institution laundries. Allows for multi-patient use. Lower wash/dry temperatures will lengthen product life. If using bleach, use lower concentrations.



Treated with AEM 5700 Antimicrobial, EPA #64881-1 ÆGIS Environmentals Midland, MI USA

The unique, patent pending constructions of DermaSaver products were developed by Ed Goodwin, founder and inventor of HipSaver hip protectors. DermaSaver products are available exclusively from HipSaver, Inc.

Call 1-800-358-4477